



Product Specification
**Microvette® CB 300 Fluoride/heparin FH, 300 µl,
cap yellow, push cap, flat base**

Page 1



Product description

Order number	16.446
Product description	Microvette® CB 300 Fluoride/heparin FH, Veterinary medicine, preparation: Fluoride + heparin, nominal volume: 300 µl, (LxØ) with cap: 46.5 x 10.8 mm, cap yellow, push cap, flat base, inner receptacle: trumpet-shaped, 100 piece(s)/bag

Product characteristics

Type of collection	capillary
Type of preparation	Fluoride + heparin
Colour code	EU
Cap	Soft stopper with retaining loop
Closure type	push cap
Application	Veterinary medicine
Base shape	flat base
Recommended centrifugation	2000 U/min
Centrifugation time	5 Min
Shape of inner tube	trumpet-shaped

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Microvette® CB 300 Fluoride/heparin FH, 300 µl,
cap yellow, push cap, flat base**

Page 2

Size

Sample volume	300 µl
Diameter	10.8 mm
Length including cap	46.5 mm
Length excluding cap	43.7 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	yellow

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	80000
Depth of case	378 mm
Width of case	288 mm
Height of case	180 mm
Case volume	0.0196 cbm
Weight of product	0.002 kg
Weight of case	2.42 kg
EAN of inner box	4038917268318
EAN case	4038917282741

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

