



Product description

Order number	17.2113.101
Product description	Minivette® POCT EDTA K3E, nominal volume: 100 µl, preparation: K3 EDTA, violet, colour code ISO, 200 piece(s)/bag, 200 piece(s)/case

Product characteristics

Type of collection	capillary
Type of preparation	K3 EDTA
Colour code	ISO
Application	capillary blood collection
Coding colour	violet

Size

Sample volume	100 µl
---------------	--------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Minivette® POCT EDTA K3E, 100 µl, plunger
violet, colour code ISO, 200 piece(s)/bag**

Page 2

Material & colours

Product material	Cellulose acetate (CA)
Colour of product	transparent

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Latex-free	yes
Batched	yes

Packaging

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	200
Piece (s) / outer case	2000
Piece(s) / pallet	144000
Depth of box	140 mm
Width of box	99 mm
Height of box	117 mm
Depth of case	513 mm
Width of case	293 mm
Height of case	140 mm
Case volume	0.021 cbm
Weight of product	0.0009 kg
Weight of case	2.51 kg
EAN of inner box	4038917375597
EAN case	4038917375580

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

