



Product description

Order number	18.1321
Product description	Microvette® CB 200 EDTA K3E, capillary blood collection, preparation: K3 EDTA, nominal volume: 200 µl, (LxØ) with cap: 46 x 10.8 mm, cap red, push cap, flat base, 100 piece(s)/bag

Product characteristics

Type of collection	capillary
Type of preparation	K3 EDTA
Colour code	EU
Closure type	push cap
Application	capillary blood collection
Base shape	flat base

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Size

Sample volume	200 µl
Diameter	10.8 mm
Length of product	40 mm
Length including cap	46 mm
Length excluding cap	40 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	red

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	64000
Depth of case	398 mm
Width of case	298 mm
Height of case	215 mm
Case volume	0.0255 cbm
Weight of product	0.0019 kg
Weight of case	2.33 kg
EAN of inner box	4038917161732
EAN case	4038917004749

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

