



Product description

Order number	18.1325.100
Product description	Microvette® CB 200 NC, capillary blood collection, preparation: Trisodium citrate solution (1:4), nominal volume: 200 µl, (LxØ) with cap: 46 x 10.8 mm, cap black, push cap, flat base, incl. sedimentation pipette, 50 piece(s)/bag

Product characteristics

Type of collection	capillary
Type of preparation	Trisodium citrate solution (1:4)
Colour code	ISO
Closure type	push cap
Application	capillary blood collection
Base shape	flat base

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Size

Sample volume	200 µl
Diameter	10.8 mm
Length of product	40 mm
Length including cap	46 mm
Length excluding cap	40 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	black

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	15000
Type of smallest subpackaging	bag
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	32000
Depth of box	271 mm
Width of box	126 mm
Height of box	36 mm
Depth of case	283 mm
Width of case	268 mm
Height of case	210 mm
Case volume	0.0159 cbm
Weight of product	0.0036 kg
Weight of case	2.6 kg
EAN of inner box	4038917127653
EAN case	4038917021333

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

