



Product description

Order number	20.1307
Product description	Microvette® 300 Fluoride/heparin FH, capillary blood collection, preparation: Fluoride + heparin, nominal volume: 300 µl, (LxØ) with cap: 47.6 x 10.8 mm, with print, cap yellow, flat base, inner receptacle: conical, 100 piece(s)/case

Product characteristics

Type of collection	capillary
Type of preparation	Fluoride + heparin
Colour code	EU
Label/ Print	with print
Colour of print/label	yellow
Cap	cap
Application	capillary blood collection
Base shape	flat base
Recommended centrifugation	2000 U/min
Centrifugation time	5 Min
Shape of inner tube	conical

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Microvette® 300 Fluoride/heparin FH, 300 µl, cap
yellow, flat base**

Page 2

Size

Sample volume	300 µl
Diameter	10.8 mm
Length including cap	47.6 mm
Length excluding cap	46.6 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	yellow

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	case
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	121000
Depth of box	139 mm
Width of box	108 mm
Height of box	66 mm
Depth of case	338 mm
Width of case	230 mm
Height of case	160 mm
Case volume	0.0124 cbm
Weight of product	0.0019 kg
Weight of case	2.43 kg
EAN of inner box	4038917086127
EAN case	4038917004701

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

