



Product Specification
Microvette® 500 Serum Gel CAT, 500 µl, cap brown, flat base

Page 1



Product description

Order number	20.1344
Product description	Microvette® 500 Serum Gel CAT, capillary blood collection, preparation: Clotting activator/gel, nominal volume: 500 µl, (LxØ) with cap: 47.6 x 10.8 mm, with print, cap brown, flat base, inner receptacle: cylindrical, 100 piece(s)/case

Product characteristics

Type of collection	capillary
Type of preparation	Clotting activator/gel
Colour code	EU/ISO
Label/ Print	with print
Colour of print/label	brown
Cap	cap
Application	capillary blood collection
Base shape	flat base
Volume of work	325 µl - 500 µl
Recommended centrifugation	10000 x g
Centrifugation time	5 Min
Shape of inner tube	cylindrical

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Microvette® 500 Serum Gel CAT, 500 µl, cap
brown, flat base**

Page 2

Size

Sample volume	500 µl
Diameter	10.8 mm
Length including cap	47.6 mm
Length excluding cap	46.6 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	brown

Purity & certification

Product category	In vitro diagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	case
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	121000
Depth of box	139 mm
Width of box	108 mm
Height of box	66 mm
Depth of case	338 mm
Width of case	228 mm
Height of case	160 mm
Case volume	0.0123 cbm
Weight of product	0.0022 kg
Weight of case	2.66 kg
EAN of inner box	4038917086264
EAN case	4038917004848

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-04

This document was prepared by EDP support and is valid without signature.



SARSTEDT