



Product description

Order number	20.1345.100
Product description	Microvette® 500 Lithium heparin LH, capillary blood collection, preparation: Lithium heparin, nominal volume: 500 µl, (LxØ) with cap: 47.6 x 10.8 mm, with print, cap green, flat base, inner receptacle: cylindrical, 100 piece(s)/case

Product characteristics

Type of collection	capillary
Type of preparation	Lithium heparin
Colour code	ISO
Label/ Print	with print
Colour of print/label	green
Cap	cap
Application	capillary blood collection
Base shape	flat base
Volume of work	325 µl - 500 µl
Recommended centrifugation	2000 x g
Centrifugation time	5 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Microvette® 500 Lithium heparin LH, 500 µl, cap green, flat base

Page 2

Shape of inner tube cylindrical

Size

Sample volume 500 µl
Diameter 10.8 mm
Length including cap 47.6 mm
Length excluding cap 46.6 mm

Material & colours

Product material Polypropylene (PP)
Colour of product transparent
Closure material High Density Polyethylene (HD-PE)
Colour of cap green

Purity & certification

Product category In vitrodiagnostic device, CE
CE certified CE
Batched yes

Packaging

Minimum order qty. 1000
Type of smallest case
subpackaging
Piece(s) / inner box 100
Piece (s) / outer case 1000
Piece(s) / pallet 121000
Depth of box 139 mm
Width of box 108 mm
Height of box 66 mm
Depth of case 338 mm
Width of case 228 mm
Height of case 160 mm
Case volume 0.0123 cbm
Weight of product 0.002 kg
Weight of case 2.49 kg
EAN of inner box 4038917127851
EAN case 4038917021432

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-08

This document was prepared by EDP support and is valid without signature.



SARSTEDT