



## Product description

Order number	20.1346.100
Product description	Microvette® 500 Lithium heparin gel LH, capillary blood collection, preparation: Lithium heparin gel, nominal volume: 500 µl, (LxØ) with cap: 47.6 x 10.8 mm, with print, cap green, flat base, inner receptacle: cylindrical, 100 piece(s)/case

## Product characteristics

Type of collection	capillary
Type of preparation	Lithium heparin gel
Colour code	ISO
Label/ Print	with print
Colour of print/label	green
Cap	cap
Application	capillary blood collection
Base shape	flat base
Volume of work	325 µl - 500 µl
Recommended centrifugation	10000 x g
Centrifugation time	5 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Microvette® 500 Lithium heparin gel LH, 500 µl,  
cap green, flat base**

**Page 2**

Shape of inner tube                      cylindrical

**Size**

Sample volume                      500 µl  
Diameter                              10.8 mm  
Length including cap                47.6 mm  
Length excluding cap                46.6 mm

**Material & colours**

Product material                      Polypropylene (PP)  
Colour of product                      transparent  
Closure material                      High Density Polyethylene (HD-PE)  
Colour of cap                              green

**Purity & certification**

Product category                      In vitrodiagnostic device, CE  
CE certified                              CE  
Batched                                  yes

**Packaging**

Minimum order qty.                      1000  
Type of smallest                              case  
subpackaging  
Piece(s) / inner box                      100  
Piece (s) / outer case                      1000  
Piece(s) / pallet                              121000  
Depth of box                              139 mm  
Width of box                              108 mm  
Height of box                              66 mm  
Depth of case                              338 mm  
Width of case                              228 mm  
Height of case                              160 mm  
Case volume                              0.0123 cbm  
Weight of product                              0.0022 kg  
Weight of case                              2.66 kg  
EAN of inner box                              4038917127868  
EAN case                                  4038917028790

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-26

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**