



Product description

Order number	20.1346
Product description	Microvette® 500 Lithium heparin gel LH, capillary blood collection, preparation: Lithium heparin gel, nominal volume: 500 µl, (LxØ) with cap: 47.6 x 10.8 mm, with print, cap orange, flat base, inner receptacle: cylindrical, 100 piece(s)/case

Product characteristics

Type of collection	capillary
Type of preparation	Lithium heparin gel
Colour code	EU
Label/ Print	with print
Colour of print/label	orange
Cap	cap
Application	capillary blood collection
Base shape	flat base
Volume of work	325 µl - 500 µl
Recommended centrifugation	10000 x g
Centrifugation time	5 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Microvette® 500 Lithium heparin gel LH, 500 µl,
cap orange, flat base**

Page 2

Shape of inner tube cylindrical

Size

Sample volume 500 µl
Diameter 10.8 mm
Length including cap 47.6 mm
Length excluding cap 46.6 mm

Material & colours

Product material Polypropylene (PP)
Colour of product transparent
Closure material High Density Polyethylene (HD-PE)
Colour of cap orange

Purity & certification

Product category In vitrodiagnostic device, CE
CE certified CE
Batched yes

Packaging

Minimum order qty. 1000
Type of smallest case
subpackaging
Piece(s) / inner box 100
Piece (s) / outer case 1000
Piece(s) / pallet 121000
Depth of box 139 mm
Width of box 108 mm
Height of box 66 mm
Depth of case 338 mm
Width of case 228 mm
Height of case 160 mm
Case volume 0.0123 cbm
Weight of product 0.0022 kg
Weight of case 2.66 kg
EAN of inner box 4038917086288
EAN case 4038917004862

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

