




**Product Specification**  
**Micro sample tube 100 µl brilliant, 100 µl, screw cap**

**Page 1**



**Product description**

Order number	42.117
Product description	Micro sample tube 100 µl brilliant, preparation: 1% brilliant cresyl blue solution, nominal volume: 100 µl, (LxØ) with cap: 46 x 10.8 mm, screw cap, cap: blue, conical base, 50 piece(s)/case

**Product characteristics**

Type of collection	Venous and capillary
Type of preparation	1% brilliant cresyl blue solution
Preparation concentration	1% brilliant cresyl blue solution
Closure type	screw cap
Example of use	Reticulocyte counting
Base shape	conical base
Volume of work	100 µl

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Micro sample tube 100 µl brilliant, 100 µl, screw cap**

**Page 2**

**Size**

Sample volume	100 µl
Diameter	10.8 mm
Length of product	46 mm
Length including cap	46 mm
Length excluding cap	44 mm

**Material & colours**

Colour of product	transparent
Closure material	Polypropylene (PP)
Colour of cap	blue

**Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

**Packaging**

Minimum order qty.	500
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	117000
Depth of case	283 mm
Width of case	173 mm
Height of case	135 mm
Case volume	0.0066 cbm
Weight of product	0.0015 kg
Weight of case	0.92 kg
EAN of inner box	4038917080897
EAN case	4038917365130

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

