



Product Specification

**Sample tube, SARSTEDT haemolysis solution,
1,000 µl, cap white, (LxØ): 97 x 15.3 mm, with print**

Page 1



Product description

Order number	45.498.003
Product description	Sample tube, Blood sugar determination from capillary blood, nominal volume: 1,000 µl, cap white, preparation: SARSTEDT haemolysis solution, screw cap, (LxØ): 97 x 15.3 mm, conical false bottom, rounded tube bottom, with print, 100 piece(s)/bag

Product characteristics

Type of collection	capillary
Type of preparation	SARSTEDT haemolysis solution
Label/ Print	with print
Colour of print/label	white
Closure type	screw cap
Application	capillary blood collection
Example of use	Blood sugar determination from capillary blood
Base shape	conical false bottom, rounded tube bottom
Volume of work	20 µl
Shape of inner tube	conical

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

**Sample tube, SARSTEDT haemolysis solution,
1,000 µl, cap white, (LxØ): 97 x 15.3 mm, with print**

Page 2

Size

Sample volume	1.000.00 µl
Diameter	15.3 mm
Length of product	97 mm
Length including cap	97 mm
Length excluding cap	92 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	white
Tube material	Polypropylene (PP)

Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In vitro diagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	500
Piece(s) / pallet	32000
Depth of case	398 mm
Width of case	298 mm
Height of case	215 mm
Case volume	0.0255 cbm
Weight of product	0.0061 kg
Weight of case	3.44 kg
EAN of inner box	4038917273022
EAN case	4038917273039

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

