



Product description

Order number	95.1048
Product description	Capillary holder, capillary blood collection, 5 piece(s)/bag

Product characteristics

Type of collection	capillary
Application	capillary blood collection

Material & colours

Colour of product	red
-------------------	-----

Purity & certification

Product category	no medical device   IVD
Batched	no

Packaging

Minimum order qty.	5
Type of smallest subpackaging	bag
Piece(s) / inner box	5
Piece (s) / outer case	1000

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

### Capillary holder

Page 2

Piece(s) / pallet	108000
Depth of case	308 mm
Width of case	223 mm
Height of case	195 mm
Case volume	0.0134 cbm
Weight of product	0.0019 kg
Weight of case	2.31 kg
EAN of inner box	4038917085120
EAN case	4038917003810

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-08

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**