



Product Specification
Faeces tube, with spoon, push cap, (LxØ): 75 x 23.5 mm, transparent

Page 1



Product description

Order number	80.620
Product description	Faeces tube, with spoon, transparent, push cap, cap: brown, (LxØ) without cap: 75 x 23.5 mm, 50 piece(s)/bag

Product characteristics

Closure type	push cap
Design	spoon

Size

Diameter	23.5 mm
Length of product	87 mm
Length excluding cap	75 mm

Material & colours

Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	brown
Tube material	Polystyrene (PS)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Faeces tube, with spoon, push cap, (LxØ): 75 x 23.5 mm, transparent

Page 2

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	18000
Depth of case	398 mm
Width of case	398 mm
Height of case	310 mm
Case volume	0.0491 cbm
Weight of product	0.0091 kg
Weight of case	5.29 kg
EAN of inner box	4038917082952
EAN case	4038917001939

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

