



Product Specification
Faeces tube, with spoon, push cap, (LxØ): 75 x 23.5 mm, transparent, sterile

Page 1



Product description

Order number	80.621
Product description	Faeces tube, with spoon, transparent, push cap, cap: brown, (LxØ) without cap: 75 x 23.5 mm, sterile, 50 piece(s)/bag

Product characteristics

Closure type	push cap
Design	spoon

Size

Diameter	23.5 mm
Length of product	87 mm
Length excluding cap	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

Faeces tube, with spoon, push cap, (LxØ): 75 x 23.5 mm, transparent, sterile

Page 2

Material & colours

Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	brown
Tube material	Polystyrene (PS)

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order qty.	250
Type of smallest subpackaging	bag
Piece(s) / inner box	50
Piece (s) / outer case	250
Piece(s) / pallet	13000
Depth of case	548 mm
Width of case	368 mm
Height of case	130 mm
Case volume	0.0262 cbm
Weight of product	0.0091 kg
Weight of case	2.73 kg
EAN of inner box	4038917082969
EAN case	4038917001946

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-05

This document was prepared by EDP support and is valid without signature.



SARSTEDT