



Product Specification

Faeces tube, with spoon, screw cap, (LxØ): 107 x 25 mm, transparent, sterile

Page 1



Product description

Order number	80.622
Product description	Faeces tube, with spoon, transparent, screw cap, Screw cap with stool spoon holder, assembled, cap: brown, (LxØ) without cap: 107 x 25 mm, sterile, 50 piece(s)/bag

Product characteristics

Cap	Screw cap with stool spoon holder
Closure type	screw cap
Design	spoon
Base shape	flat-conical base

Size

Diameter	25 mm
Length of product	111 mm
Length excluding cap	107 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

Faeces tube, with spoon, screw cap, (LxØ): 107 x 25 mm, transparent, sterile

Page 2

Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	brown
Stem material	Polystyrene (PS)
Tube material	Polypropylene (PP)

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order qty.	250
Type of smallest subpackaging	bag
Piece(s) / inner box	50
Piece (s) / outer case	250
Piece(s) / pallet	10000
Depth of case	598 mm
Width of case	398 mm
Height of case	158 mm
Case volume	0.0375 cbm
Weight of product	0.0129 kg
Weight of case	4.29 kg
EAN of inner box	4038917082976
EAN case	4066916010805

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-13

This document was prepared by EDP support and is valid without signature.



SARSTEDT