



## Product Specification

Faeces tube, with spoon, screw cap, (LxØ): 107 x 25 mm, transparent

Page 1



### Product description

Order number	80.622.111
Product description	Faeces tube, with spoon, transparent, screw cap, Screw cap with stool spoon, assembled, cap: brown, (LxØ) without cap: 107 x 25 mm, 250 piece(s)/bag

### Product characteristics

Cap	Screw cap with stool spoon
Closure type	screw cap
Design	spoon
Base shape	flat-conical base

### Size

Diameter	25 mm
Length of product	111 mm
Length excluding cap	107 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

**Faeces tube, with spoon, screw cap, (LxØ): 107 x 25 mm, transparent**

**Page 2**

### Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	brown
Stem material	Polystyrene (PS)
Tube material	Polypropylene (PP)

### Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

### Packaging

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	250
Piece (s) / outer case	500
Piece(s) / pallet	12000
Depth of case	398 mm
Width of case	398 mm
Height of case	415 mm
Case volume	0.0657 cbm
Weight of product	0.0129 kg
Weight of case	7.14 kg
EAN of inner box	4038917115346
EAN case	4038917014533

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2025-12-14

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**