



**Product Specification**  
**Faeces tube, with blade, screw cap, (LxØ): 54 x 28 mm, white**

**Page 1**



**Product description**

Order number	80.9924.014
Product description	Faeces tube, with blade, white, screw cap, cap: brown, (LxØ) without cap: 54 x 28 mm, with paper label, 50 piece(s)/bag

**Product characteristics**

Label/ Print	with paper label
Colour of print/label	black
Closure type	screw cap
Design	blade

**Size**

Diameter	28 mm
Length of product	55 mm
Length excluding cap	54 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

**Faeces tube, with blade, screw cap, (LxØ): 54 x 28 mm, white**

**Page 2**

### Material & colours

Colour of product	white
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	brown
Tube material	Polypropylene (PP)

### Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

### Packaging

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	14000
Depth of case	458 mm
Width of case	238 mm
Height of case	373 mm
Case volume	0.0407 cbm
Weight of product	0.0084 kg
Weight of case	4.9 kg
EAN of inner box	4038917154772
EAN case	4038917039598

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2025-12-15

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**