



Product description

| | |
|---------------------|---|
| Order number | 93.851.174 |
| Product description | Block Rack D17, for 50 collection tubes, Ø opening: 17 mm, format: 5 x 10, blue, with handle, material: PP |

Product characteristics

| | |
|------------------------|--------|
| Number of wells | 50 |
| Storage capacity | 50 |
| Format | 5 x 10 |
| Autoclavable (max. °C) | 121 °C |

Size

| | |
|---------------------|--------|
| Diameter of opening | 17 mm |
| Width of product | 109 mm |
| Height of product | 45 mm |
| Length of product | 246 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Block Rack D17, Ø opening: 17 mm, 5 x 10, blue,
with handle**

Page 2

Material & colours

| | |
|-------------------|--------------------|
| Product material | Polypropylene (PP) |
| Colour of product | blue |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Batched | no |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1 |
| Type of smallest subpackaging | - |
| Piece(s) / inner box | 1 |
| Piece (s) / outer case | 12 |
| Piece(s) / pallet | 840 |
| Depth of case | 453 mm |
| Width of case | 258 mm |
| Height of case | 165 mm |
| Case volume | 0.0193 cbm |
| Weight of product | 0.1983 kg |
| Weight of case | 2.79 kg |
| EAN of inner box | 4038917121866 |
| EAN case | 4038917075008 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-30

This document was prepared by EDP support and is valid without signature.



SARSTEDT