



Product Specification
**Block Rack D17, Ø opening: 17 mm, 5 x 10,
creamy white**

Page 1



Product description

Order number	93.852.170
Product description	Block Rack D17, for 50 collection tubes, Ø opening: 17 mm, format: 5 x 10, creamy white, material: PP

Product characteristics

Label/ Print	with print
Number of wells	50
Storage capacity	50
Format	5 x 10
Autoclavable (max. °C)	121 °C

Size

Diameter of opening	17 mm
Width of product	109 mm
Height of product	45 mm
Length of product	209 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2024-03-28

This document was prepared by EDP
support and is valid without signature.



SARSTEDT



Product Specification
**Block Rack D17, Ø opening: 17 mm, 5 x 10,
creamy white**

Page 2

Material & colours

Product material	Polypropylene (PP)
Colour of product	creamy white

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	1
Type of smallest subpackaging	-
Piece(s) / inner box	1
Piece (s) / outer case	20
Piece(s) / pallet	760
Depth of case	568 mm
Width of case	433 mm
Height of case	120 mm
Case volume	0.0295 cbm
Weight of product	0.1958 kg
Weight of case	4.62 kg
EAN of inner box	4038917122054
EAN case	4038917044905

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2024-03-28

This document was prepared by EDP support and is valid without signature.



SARSTEDT