



Product description

| | |
|---------------------|---|
| Order number | 93.852.171 |
| Product description | Block Rack D17, for 50 collection tubes, Ø opening: 17 mm, format: 5 x 10, red, material: recycled PP |

Product characteristics

| | |
|------------------------|------------|
| Label/ Print | with print |
| Number of wells | 50 |
| Storage capacity | 50 |
| Format | 5 x 10 |
| Autoclavable (max. °C) | 121 °C |

Size

| | |
|---------------------|--------|
| Diameter of opening | 17 mm |
| Width of product | 109 mm |
| Height of product | 45 mm |
| Length of product | 209 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

| | |
|-------------------|--------------------------------------|
| Product material | recycled Polypropylene (recycled PP) |
| Colour of product | red |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Batched | no |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1 |
| Type of smallest subpackaging | - |
| Piece(s) / inner box | 1 |
| Piece (s) / outer case | 20 |
| Piece(s) / pallet | 760 |
| Depth of case | 568 mm |
| Width of case | 433 mm |
| Height of case | 120 mm |
| Case volume | 0.0295 cbm |
| Weight of product | 0.1958 kg |
| Weight of case | 4.62 kg |
| EAN of inner box | 4038917122061 |
| EAN case | 4038917018128 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

