



## Product description

Order number	93.852.173
Product description	Block Rack D17, for 50 collection tubes, Ø opening: 17 mm, format: 5 x 10, green, material: recycled PP

## Product characteristics

Label/ Print	with print
Number of wells	50
Storage capacity	50
Format	5 x 10
Autoclavable (max. °C)	121 °C

## Size

Diameter of opening	17 mm
Width of product	109 mm
Height of product	45 mm
Length of product	209 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	recycled Polypropylene (recycled PP)
Colour of product	green

### **Purity & certification**

Product category	no medical device   IVD
Batched	no

### **Packaging**

Minimum order qty.	1
Type of smallest subpackaging	-
Piece(s) / inner box	1
Piece (s) / outer case	20
Piece(s) / pallet	760
Depth of case	568 mm
Width of case	433 mm
Height of case	120 mm
Case volume	0.0295 cbm
Weight of product	0.1958 kg
Weight of case	4.62 kg
EAN of inner box	4038917122085
EAN case	4038917018142

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

