



Product description

Order number	93.852.300
Product description	Double block rack D17, for 100 collection tubes, Ø opening: 17 mm, format: 5 x 20, grey, material: PP

Product characteristics

Number of wells	100
Storage capacity	100
Format	5 x 20
Autoclavable (max. °C)	121 °C

Size

Diameter of opening	17 mm
Width of product	109 mm
Height of product	45 mm
Length of product	418 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Double block rack D17, Ø opening: 17 mm, 5 x 20, grey

Material & colours

Product material	Polypropylene (PP)
Colour of product	grey

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	1
Type of smallest subpackaging	-
Piece(s) / inner box	1
Piece (s) / outer case	12
Piece(s) / pallet	528
Depth of case	453 mm
Width of case	433 mm
Height of case	165 mm
Case volume	0.0324 cbm
Weight of product	0.3917 kg
Weight of case	5.41 kg
EAN of inner box	4038917122221
EAN case	4038917018296

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

