



### **Product description**

Order number	93.853.300
Product description	Double block rack D13, for 100 collection tubes, Ø opening: 13 mm, format: 5 x 20, grey, material: PP

### **Product characteristics**

Number of wells	100
Storage capacity	100
Format	5 x 20
Autoclavable (max. °C)	121 °C

### **Size**

Diameter of opening	13 mm
Width of product	96 mm
Height of product	45 mm
Length of product	360 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Double block rack D13, Ø opening: 13 mm, 5 x 20, grey**

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	grey

**Purity & certification**

Product category	not regulated
CE certified	no CE mark
Batched	no

**Packaging**

Minimum order qty.	1
Type of smallest subpackaging	-
Piece(s) / inner box	1
Piece (s) / outer case	8
Piece(s) / pallet	720
Depth of case	398 mm
Width of case	368 mm
Height of case	120 mm
Case volume	0.0176 cbm
Weight of product	0.3509 kg
Weight of case	3.27 kg
EAN of inner box	4038917122634
EAN case	4038917074087

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

