



Product Specification
**Block Rack D12, Ø opening: 12 mm, 5 x 10,
orange**

Page 1



Product description

Order number	93.854.125
Product description	Block Rack D12, for 50 collection tubes, Ø opening: 12 mm, format: 5 x 10, orange, material: PP

Product characteristics

Label/ Print	with print
Number of wells	50
Storage capacity	50
Format	5 x 10
Autoclavable (max. °C)	121 °C

Size

Diameter of opening	12 mm
Width of product	80 mm
Height of product	45 mm
Length of product	152 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Block Rack D12, Ø opening: 12 mm, 5 x 10,
orange**

Page 2

Material & colours

Product material	Polypropylene (PP)
Colour of product	orange

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

Packaging

Minimum order qty.	1
Type of smallest subpackaging	-
Piece(s) / inner box	1
Piece (s) / outer case	20
Piece(s) / pallet	1560
Depth of case	430 mm
Width of case	325 mm
Height of case	130 mm
Case volume	0.0182 cbm
Weight of product	0.1214 kg
Weight of case	2.8 kg
EAN of inner box	4038917122689
EAN case	4038917018524

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-13

This document was prepared by EDP support and is valid without signature.



SARSTEDT