



Product Specification

Tube with false bottom, 5 ml, (LxØ): 92 x 15.3 mm,
conical false bottom, flat tube bottom, PP, without
cap, 100 piece(s)/bag

Page 1



Product description

Order number	60.611.010
Product description	Tube with false bottom, working volume: 5 ml, (LxØ): 92 x 15.3 mm, conical false bottom, flat tube bottom, transparent, material: PP, without cap, 100 piece(s)/bag, 1,000 piece(s)/carton

Product characteristics

Graduation	no
Base shape	conical false bottom, flat tube bottom
Volume of work	5 ml

Size

Diameter	15.3 mm
Length of product	92 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-15

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

Tube with false bottom, 5 ml, (LxØ): 92 x 15.3 mm, conical false bottom, flat tube bottom, PP, without cap, 100 piece(s)/bag

Page 2

Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Satisfies the requirement IATA, ADR
Product category In-vitro diagnostic (IVD)
CE certified CE - manufacturer's self-declaration
Batched yes

Packaging

Minimum order qty. 1000
Type of smallest subpackaging bag
Piece(s) / inner box 100
Piece (s) / outer case 1000
Piece(s) / pallet 36000
Depth of case 398 mm
Width of case 398 mm
Height of case 310 mm
Case volume 0.0491 cbm
Weight of product 0.004 kg
Weight of case 4.68 kg
EAN of inner box 4038917114820
EAN case 4038917014427

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-15

This document was prepared by EDP support and is valid without signature.



SARSTEDT