




Product Specification

**Screw cap tube, 3.5 ml, (LxØ): 92 x 13 mm,
conical false bottom, rounded tube bottom, PP,
without cap, 100 piece(s)/bag**

Page 1



Product description

Order number	60.617.010
Product description	Screw cap tube, working volume: 3.5 ml, (LxØ): 92 x 13 mm, conical false bottom, rounded tube bottom, transparent, material: PP, without cap, 100 piece(s)/bag, 1,000 piece(s)/carton

Product characteristics

Graduation	no
Base shape	conical false bottom, rounded tube bottom
Volume of work	3,5 ml

Size

Diameter	13 mm
Length of product	92 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

**Screw cap tube, 3.5 ml, (LxØ): 92 x 13 mm,
conical false bottom, rounded tube bottom, PP,
without cap, 100 piece(s)/bag**

Page 2

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	48000
Depth of case	390 mm
Width of case	377 mm
Height of case	220 mm
Case volume	0.0323 cbm
Weight of product	0.0026 kg
Weight of case	3.37 kg
EAN of inner box	4038917115094
EAN case	4038917079167

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-15

This document was prepared by EDP support and is valid without signature.



SARSTEDT