



Product description

| | |
|---------------------|---|
| Order number | 75.1354.001 |
| Product description | Container with screw cap, collection and storage of urine, max. working volume: 100 ml, Ø opening: 58 mm, transparent, graduated, material: PP, screw cap, cap: green, cap enclosed, 25 pieces each stacked inside each other, 500 piece(s)/bag |

Product characteristics

| | |
|----------------|---------------------------------|
| Graduation | yes |
| Closure type | screw cap |
| Application | collection and storage of urine |
| Volume of work | 100 ml |

Size

| | |
|----------------------|--------|
| Sample volume | 100 ml |
| Max. volume | 100 ml |
| Diameter of opening | 58 mm |
| Length including cap | 73 mm |
| Length excluding cap | 71 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

| | |
|-------------------|-----------------------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |
| Closure material | High Density Polyethylene (HD-PE) |
| Colour of cap | green |

Purity & certification

| | |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified | CE supplier |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 500 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 500 |
| Piece (s) / outer case | 500 |
| Piece(s) / pallet | 12000 |
| Depth of case | 600 mm |
| Width of case | 330 mm |
| Height of case | 400 mm |
| Case volume | 0.0778 cbm |
| Weight of product | 0.0115 kg |
| Weight of case | 7.64 kg |
| EAN of inner box | 4038917247801 |
| EAN case | 8052109529632 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

