



Product Specification
**Urine container, 75 ml, (ØxH): 65 x 48 mm, PP,
transparent**

Page 1



Product description

Order number	75.1437.100
Product description	Urine container, max. working volume: 75 ml, (ØxH): 65 x 48 mm, transparent, graduated, material: PP

Product characteristics

Graduation	yes
Volume of work	75 ml

Size

Sample volume	75 ml
Max. volume	75 ml
Diameter	65 mm
Height of product	48 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Urine container, 75 ml, (ØxH): 65 x 48 mm, PP,
transparent**

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	1000
Piece(s) / inner box	25
Piece (s) / outer case	1000
Piece(s) / pallet	28000
Case volume	0.065 cbm
Weight of product	0.0092 kg
Weight of case	9.65 kg
EAN of inner box	4038917129138
EAN case	4038917074711

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

