



Product description

Order number	75.560.901
Product description	Urine container, max. working volume: 148 ml, (ØxH): 73 x 67 mm, white, material: PS

Product characteristics

Graduation	no
Volume of work	148 ml

Size

Sample volume	148 ml
Diameter	73 mm
Height of product	67 mm

Material & colours

Product material	Polystyrene (PS)
Colour of product	white

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	2500
Piece(s) / inner box	100
Piece (s) / outer case	2500
Piece(s) / pallet	90000
Depth of case	457 mm
Width of case	368 mm
Height of case	356 mm
Case volume	0.0599 cbm
Weight of product	0.0028 kg
Weight of case	8.16 kg
EAN of inner box	4038917318198
EAN case	4038917916547

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

