



Product Specification

Container with screw cap, 100 ml, (ØxH): 57 x 76 mm, PP, with safety label, transparent

Page 1



Product description

Order number

75.562.105

Product description

Container with screw cap, collection and storage of urine, max. working volume: 100 ml, (ØxH): 57 x 76 mm, Ø opening: 62 mm, transparent, graduated, material: PP, screw cap, cap: white, with safety label, cap assembled, sterile, 5 piece(s)/bag

Product characteristics

Label/ Print

with safety label

Colour of print/label

white/red

Graduation

yes

Closure type

screw cap

Application

collection and storage of urine

Volume of work

100 ml

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Container with screw cap, 100 ml, (ØxH): 57 x 76 mm, PP, with safety label, transparent

Page 2

Size

Sample volume	100 ml
Max. volume	100 ml
Diameter	57 mm
Diameter of opening	62 mm
Height of product	76 mm
Length excluding cap	72 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	white

Purity & certification

Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order qty.	200
Type of smallest subpackaging	bag
Piece(s) / inner box	5
Piece (s) / outer case	200
Piece(s) / pallet	5600
Depth of case	588 mm
Width of case	398 mm
Height of case	255 mm
Case volume	0.0597 cbm
Weight of product	0.0155 kg
Weight of case	3.78 kg
EAN of inner box	4038917113632
EAN case	4038917014052

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

