



## Product Specification

**Multi-purpose container, max. volume: 125 ml,  
(LxØ): 85 x 62 mm, graduated, PP, transparent**

**Page 1**



### Product description

|                     |  |
|---------------------|--|
| Order number        | 75.565   |
| Product description | Multi-purpose container, max. volume: 125 ml, (LxØ): 85 x 62 mm, Ø opening: 62 mm, transparent, graduated, material: PP, 50 piece(s)/bag |

### Product characteristics

|                |                                 |
|----------------|---------------------------------|
| Graduation     | yes                             |
| Application    | collection and storage of urine |
| Volume of work | 125 ml                          |

### Size

|                      |        |
|----------------------|--------|
| Sample volume        | 125 ml |
| Max. volume          | 125 ml |
| Diameter             | 62 mm  |
| Diameter of opening  | 62 mm  |
| Height of product    | 85 mm  |
| Length of product    | 85 mm  |
| Length excluding cap | 85 mm  |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

**Multi-purpose container, max. volume: 125 ml,  
(LxØ): 85 x 62 mm, graduated, PP, transparent**

**Page 2**

### Material & colours

|                   |                    |
|-------------------|--------------------|
| Product material  | Polypropylene (PP) |
| Colour of product | transparent        |

### Purity & certification

|                  |               |
|------------------|---------------|
| Product category | not regulated |
| CE certified     | no CE mark    |
| Batched          | yes           |

### Packaging

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 500           |
| Type of smallest subpackaging | bag           |
| Piece(s) / inner box          | 50            |
| Piece (s) / outer case        | 500           |
| Piece(s) / pallet             | 18000         |
| Depth of box                  | 320 mm        |
| Width of box                  | 170 mm        |
| Height of box                 | 65 mm         |
| Depth of case                 | 578 mm        |
| Width of case                 | 368 mm        |
| Height of case                | 195 mm        |
| Case volume                   | 0.0415 cbm    |
| Weight of product             | 0.0107 kg     |
| Weight of case                | 5.93 kg       |
| EAN of inner box              | 4038917082525 |
| EAN case                      | 4038917001601 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2025-12-14

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**