



### **Product description**

Order number	76.1340.560
Product description	Lid, for 30-ml medicine cup, material: PE, white, collection and storage of urine, 75 piece(s)/bag

### **Product characteristics**

Application	collection and storage of urine
Example of use	cap for urine collection tube
Suitable for	30-ml medicine cup

### **Material & colours**

Product material	Polyethylene (PE)
Colour of product	white
Colour of cap	white

### **Purity & certification**

Product category	Medical device (MD)
CE certified	CE supplier
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## **Packaging**

Minimum order qty.	3750
Type of smallest subpackaging	bag
Piece(s) / inner box	3750
Piece (s) / outer case	3750
Piece(s) / pallet	90000
Depth of box	390 mm
Width of box	265 mm
Height of box	590 mm
Depth of case	600 mm
Width of case	395 mm
Height of case	275 mm
Case volume	0.0652 cbm
Weight of product	0.0007 kg
Weight of case	3.64 kg
EAN of inner box	4038917127806
EAN case	4042836007156

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

