



Product description

Order number	76.561
Product description	Snap-on lid, suitable for 250-ml container 75.560, material: PS, highly transparent, 50 piece(s)/tube bag

Product characteristics

Closure type	snap-on lid
Suitable for	250-ml container 75.560

Material & colours

Product material	Polystyrene (PS)
Colour of product	highly transparent

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Snap-on lid, suitable for 250-ml container 75.560,
PS, highly transparent

Page 2

Packaging

Minimum order qty.	500
Type of smallest subpackaging	tube bag
Piece(s) / inner box	500
Piece (s) / outer case	500
Piece(s) / pallet	30000
Depth of box	210 mm
Width of box	80 mm
Height of box	80 mm
Depth of case	423 mm
Width of case	213 mm
Height of case	185 mm
Case volume	0.0167 cbm
Weight of product	0.0013 kg
Weight of case	0.91 kg
EAN of inner box	4038917339926
EAN case	4038917339919

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-14

This document was prepared by EDP support and is valid without signature.



SARSTEDT