



**Product Specification**  
**S-Monovette® , 9 ml, for Formalin system, (LxØ):**  
**92 x 16 mm, with paper label**

**Page 1**



**Product description**

Order number	09.1704.001
Product description	S-Monovette® , up to 9 ml, for Formalin system, (LxØ): 92 x 16 mm, with paper label, 50 piece(s)/bag

**Product characteristics**

Type of preparation	without preparation
Label/ Print	with paper label
Colour of print/label	black
Graduation	fill mark
Closure type	screw cap
Application	Tissue fixation
Suitable for	Formalin system
Volume of work	9 ml
Liquid-tight	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**S-Monovette® , 9 ml, for Formalin system, (LxØ):**  
**92 x 16 mm, with paper label**

**Page 2**

**Size**

Max. volume	9 ml
Outer diameter of cap	19 mm
Diameter	16 mm
Length of product	92 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	white
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

**Purity & certification**

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	non-sterile
Batched	yes

**Packaging**

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	24000
Depth of case	388 mm
Width of case	378 mm
Height of case	208 mm
Case volume	0.0305 cbm
Weight of product	0.0078 kg
Weight of case	4.39 kg
EAN of inner box	4038917335317
EAN case	4038917335300

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

