



**Product Specification**  
**Formalin system, buffered formalin solution 450 ml**

**Page 1**



**Product description**

Order number	51.1703
Product description	Formalin system, preparation: buffered formalin solution 450 ml, 1 piece(s)/case

**Product characteristics**

Type of preparation	buffered formalin solution 450 ml
Label/ Print	with paper label
Application	Tissue fixation

**Size**

Max. volume	450 ml
Width of product	98 mm
Length of product	330 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Purity & certification**

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

**Packaging**

Minimum order qty.	6
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	6
Piece(s) / pallet	360
Depth of box	96 mm
Width of box	63 mm
Height of box	255 mm
Depth of case	363 mm
Width of case	325 mm
Height of case	164 mm
Case volume	0.0193 cbm
Weight of product	0.5302 kg
Weight of case	4.22 kg
EAN of inner box	4038917336031
EAN case	4038917336024

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

