



Product description

| | |
|---------------------|---|
| Order number | 51.1703.025 |
| Product description | Formalin Set, preparation: buffered formalin solution 450 ml, set components: Formalin system with 450 ml formaldehyde solution (4%) and 50 S-Monovettes 25 ml, 1 piece(s)/case |

Product characteristics

| | |
|---------------------|-----------------------------------|
| Type of preparation | buffered formalin solution 450 ml |
| Label/ Print | with paper label |
| Design | S-Monovette® 25 ml |
| Application | Tissue fixation |

Size

| | |
|-------------------|--------|
| Width of product | 98 mm |
| Length of product | 330 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Formalin Set, buffered formalin solution 450 ml,
S-Monovette® 25 ml**

Page 2

Purity & certification

| | |
|------------------|--------------------------------------|
| Product category | In-vitro diagnostic (IVD) |
| CE certified | CE - manufacturer's self-declaration |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1 |
| Type of smallest subpackaging | case |
| Piece(s) / inner box | 1 |
| Piece (s) / outer case | 1 |
| Piece(s) / pallet | 126 |
| Depth of box | 344 mm |
| Width of box | 243 mm |
| Height of box | 103 mm |
| Depth of case | 355 mm |
| Width of case | 255 mm |
| Height of case | 123 mm |
| Case volume | 0.0111 cbm |
| Weight of product | 1.232 kg |
| Weight of case | 1.69 kg |
| EAN of inner box | 4038917336093 |
| EAN case | 4038917336086 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

