



Product Specification
Gallipot, 250 ml, (LxØ): 83 x 75 mm, PP,
translucent

Page 1



Product description

Order number	75.9924.518
Product description	Gallipot, max. working volume: 250 ml, (LxØ): 83 x 75 mm, translucent, material: PP

Product characteristics

Application	Pathology storage container
Volume of work	250 ml
Autoclavable (max. °C)	121 °C

Size

Max. volume	250 ml
Diameter	75 mm
Width of product	77.2 mm
Height of product	83 mm
Length of product	83 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Gallipot, 250 ml, (LxØ): 83 x 75 mm, PP,
translucent

Page 2

Material & colours

Product material	Polypropylene (PP)
Colour of product	translucent

Purity & certification

Latex-free	yes
Pyrogen-free	no
Safety data sheet	no

Packaging

Minimum order qty.	500
Piece(s) / inner box	500
Piece (s) / outer case	500
Piece(s) / pallet	16000
Depth of case	428 mm
Width of case	372 mm
Height of case	438 mm
Case volume	0.0697 cbm
Weight of case	8.7 kg
EAN of inner box	4038917155311
EAN case	4038917039352

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2024-03-29

This document was prepared by EDP support and is valid without signature.



SARSTEDT