



Product description

Order number	95.1706
Product description	Wall holder, material: stainless steel, for Formalin system, 1 piece(s)/case

Product characteristics

Application	Tissue fixation
Suitable for	Formalin system

Size

Width of product	101 mm
Length of product	200 mm

Material & colours

Product material	stainless steel
------------------	-----------------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	yes

Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	324
Depth of case	268 mm
Width of case	191 mm
Height of case	100 mm
Case volume	0.0051 cbm
Weight of product	0.269 kg
Weight of case	0.42 kg
EAN of inner box	4038917335461
EAN case	4038917335454

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

