



### **Product description**

Order number	95.1706
Product description	Wall holder, material: stainless steel, for Formalin system, 1 piece(s)/case

### **Product characteristics**

Application	Tissue fixation
Suitable for	Formalin system

### **Size**

Width of product	101 mm
Length of product	200 mm

### **Material & colours**

Product material	stainless steel
------------------	-----------------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category

no medical device | IVD

Batched

yes

Packaging

Minimum order qty.

1

Type of smallest subpackaging

case

Piece(s) / inner box

1

Piece (s) / outer case

1

Piece(s) / pallet

324

Depth of case

268 mm

Width of case

191 mm

Height of case

100 mm

Case volume

0.0051 cbm

Weight of product

0.269 kg

Weight of case

0.42 kg

EAN of inner box

4038917335461

EAN case

4038917335454

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

