



**Product Specification**  
**Lid for sample tube 73.641 & 73.646, length: 8 mm, LD-PE**

**Page 1**



**Product description**

Order number	65.649
Product description	Lid for sample tube 73.641 & 73.646, natural, length: 8 mm, material: LD-PE, 1,000 piece(s)/bag

**Product characteristics**

Cap	Skirted push cap
Suitable for	sample tube 73.641 & 73.646

**Size**

Length of product	8 mm
-------------------	------

**Material & colours**

Product material	Low Density Polyethylene (LD-PE)
Colour of product	natural

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Lid for sample tube 73.641 & 73.646, length: 8 mm, LD-PE**

**Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

**Packaging**

Minimum order qty.	10000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	10000
Piece(s) / pallet	420000
Depth of case	388 mm
Width of case	378 mm
Height of case	250 mm
Case volume	0.0367 cbm
Weight of product	0.0005 kg
Weight of case	5.53 kg
EAN of inner box	4038917083119
EAN case	4038917002073

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

