



Product Specification
Insert, suitable for tubes and S-Monovette® Ø 13 mm, transparent

Page 1



Product description

Order number	73.1051
Product description	Insert, suitable for tubes and S-Monovette® Ø 13 mm, material: PS, 1,000 µl, length: 23.5 mm, transparent, 500 piece(s)/bag

Product characteristics

Graduation	no
Suitable for	tubes and S-Monovette® Ø 13 mm
Base shape	skirted conical base

Size

Max. volume	1000 µl
Diameter of opening	9 mm
Length of product	23.5 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Insert, suitable for tubes and S-Monovette® Ø 13
mm, transparent

Page 2

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

Packaging

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	2000
Piece(s) / pallet	156000
Depth of case	430 mm
Width of case	325 mm
Height of case	130 mm
Case volume	0.0182 cbm
Weight of product	0.0006 kg
Weight of case	1.6 kg
EAN of inner box	4038917085144
EAN case	4038917003827

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-14

This document was prepared by EDP support and is valid without signature.



SARSTEDT