



**Product Specification**  
**Insert, suitable for tubes and S-Monovette® Ø 15 mm, transparent**

**Page 1**



**Product description**

Order number	73.1052.001
Product description	Insert, suitable for tubes and S-Monovette® Ø 15 mm, material: PS, 2.5 ml, length: 32 mm, transparent, 500 piece(s)/bag

**Product characteristics**

Graduation	no
Suitable for	tubes and S-Monovette® Ø 15 mm
Base shape	skirted conical base

**Size**

Max. volume	2.5 ml
Diameter of opening	11 mm
Length of product	32 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Insert, suitable for tubes and S-Monovette® Ø 15**  
**mm, transparent**

**Page 2**

**Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent

**Purity & certification**

Product category	not regulated
CE certified	no CE mark
Batched	no

**Packaging**

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	2000
Piece(s) / pallet	156000
Depth of case	430 mm
Width of case	325 mm
Height of case	130 mm
Case volume	0.0182 cbm
Weight of product	0.0008 kg
Weight of case	1.98 kg
EAN of inner box	4038917124874
EAN case	4038917019507

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2025-12-13

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**