



Product Specification
Insert, suitable for tubes and S-Monovette® Ø 15 mm, transparent

Page 1



Product description

Order number	73.1052.001
Product description	Insert, suitable for tubes and S-Monovette® Ø 15 mm, material: PS, length: 32 mm, transparent, 500 piece(s)/bag

Product characteristics

Graduation	no
Suitable for	tubes and S-Monovette® Ø 15 mm
Base shape	skirted conical base

Size

Diameter of opening	11 mm
Length of product	32 mm

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Insert, suitable for tubes and S-Monovette® Ø 15
mm, transparent

Page 2

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 2000
Type of smallest bag
subpackaging
Piece(s) / inner box 500
Piece (s) / outer case 2000
Piece(s) / pallet 156000
Depth of case 428 mm
Width of case 323 mm
Height of case 130 mm
Case volume 0.018 cbm
Weight of product 0.0008 kg
Weight of case 1.98 kg
EAN of inner box 4038917124874
EAN case 4038917019507

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-07

This document was prepared by EDP support and is valid without signature.



SARSTEDT