



**Product Specification**  
**RIA tube, suitable for gamma counters,**  
**transparent**

**Page 1**



**Product description**

Order number	73.1055
Product description	RIA tube, suitable for gamma counters, material: PP, 0.5 ml, (LxØ): 22 x 8 mm, transparent, 1,000 piece(s)/bag

**Product characteristics**

Graduation	no
Closure type	push cap
Suitable for	gamma counters
Base shape	conical base

**Size**

Max. volume	0.5 ml
Diameter	8 mm
Diameter of opening	29 mm
Height of product	22 mm
Length of product	22 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification**  
**RIA tube, suitable for gamma counters,**  
**transparent**

**Page 2**

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

**Purity & certification**

Product category	no medical device   IVD
Batched	no

**Packaging**

Minimum order qty.	10000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	10000
Piece(s) / pallet	480000
Depth of case	388 mm
Width of case	378 mm
Height of case	215 mm
Case volume	0.0315 cbm
Weight of product	0.0005 kg
Weight of case	5.05 kg
EAN of inner box	4038917085175
EAN case	4038917003865

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2024-03-29

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**