



**Product Specification**  
**Sample tube, suitable for auto-analyzer,  
transparent**

**Page 1**



**Product description**

Order number	73.646
Product description	Sample tube, suitable for auto-analyzer, material: PS, 3.5 ml, (LxØ): 38 x 14 mm, transparent, 1,000 piece(s)/bag

**Product characteristics**

Graduation	no
Closure type	push cap
Suitable for	auto-analyzer
Base shape	conical base

**Size**

Max. volume	3.5 ml
Diameter	14 mm
Diameter of opening	8 mm
Height of product	38 mm
Length of product	38 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Sample tube, suitable for auto-analyzer,**  
**transparent**

**Page 2**

**Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent

**Purity & certification**

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

**Packaging**

Minimum order qty.	5000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	5000
Piece(s) / pallet	120000
Depth of case	398 mm
Width of case	398 mm
Height of case	415 mm
Case volume	0.0657 cbm
Weight of product	0.0018 kg
Weight of case	9.62 kg
EAN of inner box	4038917083096
EAN case	4038917002066

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2025-12-14

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**