



Product description

| | |
|---------------------|--|
| Order number | 73.650.500 |
| Product description | Sample tube, suitable for Coulter Counter, cap enclosed, material: PS, (LxØ): 60 x 34 mm, 250 piece(s)/bag |

Product characteristics

| | |
|--------------|-----------------|
| Graduation | no |
| Closure type | push cap |
| Suitable for | Coulter Counter |
| Base shape | flat base |

Size

| | |
|---------------------|-------|
| Diameter | 34 mm |
| Diameter of opening | 12 mm |
| Length of product | 60 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

| | |
|------------------|----------------------------------|
| Product material | Polystyrene (PS) |
| Closure material | Low Density Polyethylene (LD-PE) |

Purity & certification

| | |
|------------------|--------------------------------------|
| Product category | In-vitro diagnostic (IVD) |
| CE certified | CE - manufacturer's self-declaration |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 500 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 250 |
| Piece (s) / outer case | 500 |
| Piece(s) / pallet | 18000 |
| Depth of case | 398 mm |
| Width of case | 398 mm |
| Height of case | 310 mm |
| Case volume | 0.0491 cbm |
| Weight of product | 0.0077 kg |
| Weight of case | 4.51 kg |
| EAN of inner box | 4038917157841 |
| EAN case | 4038917157834 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

