



**Product Specification**  
**Sample tube, suitable for Boehringer ES 300® and ES 600®, transparent**

**Page 1**



**Product description**

Order number	73.663
Product description	Sample tube, suitable for Boehringer ES 300® and ES 600®, material: PS, 2 ml, (LxØ): 40 x 12 mm, transparent, 500 piece(s)/bag

**Product characteristics**

Graduation	no
Closure type	push cap
Suitable for	Boehringer ES 300® and ES 600®
Base shape	conical base

**Size**

Max. volume	2 ml
Diameter	12 mm
Diameter of opening	10 mm
Height of product	40 mm
Length of product	40 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Sample tube, suitable for Boehringer ES 300® and ES 600®, transparent**

**Page 2**

**Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent

**Purity & certification**

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

**Packaging**

Minimum order qty.	5000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	5000
Piece(s) / pallet	180000
Depth of case	398 mm
Width of case	398 mm
Height of case	310 mm
Case volume	0.0491 cbm
Weight of product	0.0012 kg
Weight of case	6.72 kg
EAN of inner box	4038917083195
EAN case	4038917002158

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2025-12-14

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**