



Product Specification
Sample tube, suitable for Hitachi analyser,
transparent

Page 1



Product description

Order number	73.666
Product description	Sample tube, suitable for Hitachi analyser, material: PS, 2 ml, (LxØ): 38 x 17 mm, transparent, 5,000 piece(s)/case

Product characteristics

Graduation	no
Closure type	push cap
Suitable for	Hitachi analyser
Base shape	conical base

Size

Max. volume	2 ml
Diameter	17 mm
Diameter of opening	10 mm
Height of product	38 mm
Length of product	38 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2024-03-29

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
Sample tube, suitable for Hitachi analyser,
transparent

Page 2

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	5000
Type of smallest subpackaging	case
Piece(s) / inner box	1000
Piece (s) / outer case	5000
Piece(s) / pallet	150000
Depth of case	398 mm
Width of case	398 mm
Height of case	340 mm
Case volume	0.0539 cbm
Weight of product	0.0017 kg
Weight of case	9.3 kg
EAN of inner box	4038917083225
EAN case	4038917002172

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2024-03-29

This document was prepared by EDP support and is valid without signature.



SARSTEDT