




Product Specification
Haemo-Diff® , with smear edge, for producing blood smears

Page 1



Product description

Order number	14.1218
Product description	Haemo-Diff® , with smear edge, for producing blood smears, 500 piece(s)/case

Product characteristics

Design	with smear edge
Application	Blood smear
Connection	S-Monovette®

Size

Length of product	21 mm
-------------------	-------

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	2000
Type of smallest subpackaging	case
Piece(s) / inner box	500
Piece (s) / outer case	2000
Piece(s) / pallet	216000
Depth of box	191 mm
Width of box	129 mm
Height of box	129 mm
Depth of case	398 mm
Width of case	268 mm
Height of case	150 mm
Case volume	0.016 cbm
Weight of product	0.0013 kg
Weight of case	3.27 kg
EAN of inner box	4038917085830
EAN case	4038917004466

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

