



## Product description

Order number	60.9924.952
Product description	Screw cap tube, working volume: 12 ml, (LxØ): 99 x 16 mm, material: SAN, round base, highly transparent, screw cap, yellow, cap assembled, 4 printed boxes, with graduations, sterile, 4 piece(s)/blister

## Product characteristics

Graduation	yes
Closure type	screw cap
Base shape	round base
Volume of work	12 ml

## Size

Max. volume	12 ml
Diameter	16 mm
Diameter of opening	14 mm
Height of product	99 mm
Length of product	99 mm
Length excluding cap	99 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Material & colours**

Product material	Styrene acrylonitrile (SAN)
Colour of product	highly transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	yellow

### **Purity & certification**

Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

### **Packaging**

Minimum order qty.	250
Type of smallest subpackaging	blister
Piece(s) / inner box	250
Piece (s) / outer case	250
Piece(s) / pallet	4000
Depth of case	548 mm
Width of case	408 mm
Height of case	373 mm
Case volume	0.0834 cbm
Weight of product	0.0257 kg
Weight of case	7.2 kg
EAN of inner box	4038917155472
EAN case	4038917357937

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

